


SYRINGE PUMP CONTROL SOFTWARE & CABLES



SyringePumpPro
Control Software for Syringe Pumps

Very Important: When purchasing SyringePumpPro please provide us with the email address who will be downloading the software.

Pumps per SyringePumpPro Software License	Model #
1	SYRINGE-PUMP-PRO-01
2	SYRINGE-PUMP-PRO-02
3	SYRINGE-PUMP-PRO-03
5	SYRINGE-PUMP-PRO-05
6+	SYRINGE-PUMP-PRO-06

Control programmable infusion syringe pumps including Microfluidic, Peristaltic and OEM models.

SyringePumpPro is an easy to use computer interface for programmable syringe pumps. It permits faster programming and operation of one or more pumps. You can create pump programs and upload them to your pumps.

Experience friendly effective software support like never before.

- Synchronously start and stop multiple pumps
- Configure and control multi-pump protocols
- Eliminates hours of manual pump operation
- Monitor pump parameters
- Upload pump programs
- End front panel button mashing
- OEM Pump? No Keypad?
- Send commands to one or many pumps
- No more re-typing of commands in terminal software
- Quick and easy to use interface for communicating with your pumps

Compatible with following brands



RS-232 TO USB CONVERTER



DESCRIPTION:

- Adds an RS-232 serial port to your computer by attaching to a USB port
- Includes software driver CD
- Attaches to cable CBL-PC-PUMP-7



PUMP-TO-PC RS-232 PRIMARY NETWORK CABLE

DESCRIPTION:

- Allows full control of pump or other device from a computer
- DB-9 adapter connects to 9-pin serial port (25-pin adapters available upon request)
- Utilizes the RS-232 communications port of device
- Available in 7' and 25' length cable
- Quick and easy setup



MODELS:

- CBL-PC-PUMP-7 (7 FT. CABLE)
- CBL-PC-PUMP-25 (25 FT. CABLE)

KFTechnology
LIFE SCIENCE TOOLS

KF Technology Srl - Via Amedeo Bocchi, 84/8 - 00125 Roma - Italy
Ph: +39 06.454.34.179 - Mobile: +39 339.533.03.22 - Fax +39 06.972.53.131 - Skype: kft2002
info@kftechnology.it - www.kftechnology.it